

Product Specifications

CL07.0227

PHARMEX

Buffer solution pH 2.00

Puffer Lösung pH 2.00
Tampon solution pH 2.00
Tampone soluzione pH 2.00
Tampón solución pH 2.00
Buffer oplossing pH 2.00

Reagent acc. Ph. Eur. 9(4.1.3) - (4000200)
6.57 g KCl + 119 ml HCl 0.1 mol/l / l H₂O

Mol weight:
Density: 1,01 g/ml
HS Nr: 38229000
CAS Nr:
UN:
ADR:

Specifications: