

Product Specifications

CL00.0121

ANALYT

Aluminium(III) nitrate.9aq a.r.

Aluminio(III) nitrato.9aq p.a.
Alluminio(III) nitrato.9aq p.a.
Aluminium(III)nitrat.9aq z.A.
Aluminium(III) nitrate.9aq p.a.
Aluminium(III)nitraat.9aq p.a.

For laboratory use, ACS
98.5+% Al(NO₃)₃.9H₂O

Mol weight: 375.13 g/mol
Density: 1,00 g/cm³
HS Nr: 28342980
CAS Nr: 7784-27-2
UN: 1438
ADR: 5.1,III



Specifications:

Assay	>98.5%
Chloride (Cl)	<0.001%
Sulfate (SO ₄)	<0.005%
Lead (Pb)	<0.001%
Iron (Fe)	<0.002%
Calcium (Ca)	<0.005%
Sodium (Na)	<0.005%
Potassium (K)	<0.002%
Copper (Cu)	<0.001%
Cadmium (Cd)	<0.001%
Ammonium (NH ₄)	<0.02%
Insoluble Matter (Non Solubles)	<0.005%
Magnesium (Mg)	<0.001%
Heavy Metals as Lead (Pb)	<0.001%