

Product Specifications

CL00.0310

ANALYT

Hydrochloric acid 37% a.r.

Salzsäure 37% z.A.
Chloorwaterstofzuur 37% p.a.
Acide chlorhydrique 37% p.a.
Acido clorhídrico 37% p.a.
Acido clorídrico 37% p.a.
Acido clorídrico 37% p.a.
For laboratory use
37+% HCl

Mol weight: 36.46 g/mol
Density: 1,19 g/ml
HS Nr: 28061000
CAS Nr: 7647-01-0
UN: 1789
ADR: 8,II

**Specifications:**

Assay	>37%
Residue after Ignition	<0.0005%
Arsenic (As)	<0.000005%
Cadmium (Cd)	<0.000001%
Copper (Cu)	<0.000005%
Iron (Fe)	<0.00005%
Lead (Pb)	<0.000005%
Zinc (Zn)	<0.0001%
Ammonium (NH ₄)	<0.0003%
Free Chlorine	<0.00005%
Sulfite (SO ₃)	<0.0002%
Sulfate (SO ₄)	<0.0001%