

# Product Specifications

CL00.0470

CHROMALYT

## Dichloromethane, HPLC grade

Dichloormethaan, HPLC grade  
Diclorometano, HPLC grade  
Dichlorométhane, HPLC grade  
Diclorometano, HPLC grade  
Dichlormethan, HPLC grade

**For laboratory use, HPLC**  
**99.8+% CH<sub>2</sub>Cl<sub>2</sub> (Stabilised with Ethanol)**

Mol weight: 84.93 g/mol  
Density: 1,32 g/ml  
HS Nr: 29031200  
CAS Nr: 75-09-2  
UN: 1593  
ADR: 6.1,III



### Specifications:

Assay	>99.8%
Water	<100ppm
Acidity	<0.0005meq/g
Non Volatiles	<0.0005%
Transmittance @ 235nm	> 40%
Transmittance @ 240nm	> 75%
Transmittance @ 250nm	> 98%
Transmittance @ 260nm	> 99%
Residue after Ignition	< 2 ppm
Filtered 0,22µm (Micron)	