

Product Specifications

CL00.0817

CHROMALYT

Hexane-(n) 99+%, HPLC grade

Hexano-(n), HPLC grade
Esano-(n), HPLC grade
Hexane-(n), HPLC grade
Hexan-(n), HPLC grade
Hexaan-(n), HPLC grade
Hesano-(n), HPLC grade
**For laboratory use, HPLC
99+% C₆H₁₄**

Mol weight: 86.18 g/mol
Density: 0,66 g/ml
HS Nr: 29011000
CAS Nr: 110-54-3
UN: 1208
ADR: 3,II

**Specifications:**

Assay	>99%
Water	<100ppm
Acidity	<0.0005meq/g
Non Volatiles	<0.0005%
Transmittance @ 200nm	> 50%
Transmittance @ 210nm	> 70%
Transmittance @ 220nm	> 90%
Transmittance @ 230nm	> 98%
Transmittance @ 240nm	> 99%
Filtered 0,22µm (Micron)	