

Product Specifications

CL00.1104

Potassium hexacyanoferrate(II).3aq v.p.

Kaliumhexacyanoferrat(II).3aq Rst.
Potassio hexacianoferrato(II).3aq t.p.
Kaliumhexacyanoferraat(II).3aq z.z.
Potassium hexacyanoferrate(II).3aq t.p.
Potassio esacianoferrato(II).3aq e.p.

For laboratory use
99+% K₄(Fe(CN)₆).3H₂O

Mol weight: 422.39 g/mol
Density: 1,85 g/cm³
HS Nr: 28372000
CAS Nr: 14459-95-1
UN:
ADR:

Specifications:

Assay	>99%
Chloride (Cl)	<0.05%
Sulfate (SO ₄)	<0.05%