



CL00.1112

Potassium nitrate v.p.

Potassium nitrate t.p.
Kaliumnitrat Rst.
Potassio nitrato e.p.
Potasio nitrato t.p.
Kaliumnitraat z.z.

**For laboratory use, BP, FCC
99-100.5% KNO₃**

Mol weight: 101.10 g/mol
Density: 2,11 g/cm³
HS Nr: 28342100
CAS Nr: 7757-79-1
UN: 1486
ADR: 5.1,III



Specifications:

Assay	: 99-100.5%
Loss on drying	<0.5%
Ammonium (NH ₄)	<0.005%
Arsenic (As)	<0.0001%
Copper (Cu)	<0.001%
Iron (Fe)	<0.001%
Sodium (Na)	<0.1%
Lead (Pb)	<0.001%
Zinc (Zn)	<0.001%
Chloride (Cl)	<0.001%
Nitrite (NO ₂)	<0.001%
Sulfate (SO ₄)	<0.01%