

# Product Specifications

**CL00.1113**

## Potassium permanganate v.p.

Potassio permanganato e.p.  
Kaliumpermanganat Rst.  
Potassium permanganate t.p.  
Potasio permanganato t.p.  
Kaliumpermanganaat z.z.

**For laboratory use, DAB, Ph. Eur., BP, Ph. Franç., USP  
99-100.5%  $\text{KMnO}_4$**

Mol weight: 158.03 g/mol  
Density: 2.70 g/cm<sup>3</sup>  
HS Nr: 28416100  
CAS Nr: 7722-64-7  
UN: 1490  
ADR: 5.1,II



### Specifications:

|                            |             |
|----------------------------|-------------|
| Assay                      | : 99-100.5% |
| Loss on drying             | <0.2%       |
| Chloride (Cl)              | <0.01%      |
| Sulfate (SO <sub>4</sub> ) | <0.03%      |