

Product Specifications

CL00.1143**Potassium chromate v.p.**

Potasio cromato t.p.
Potassium chromate t.p.
Kaliumchromat Rst.
Kaliumchromaat z.z.
Potassio cromato e.p.

For laboratory use
99+% K₂CrO₄

Mol weight: 194.20 g/mol
Density: 2.73 g/cm³
HS Nr: 28415000
CAS Nr: 7789-00-6
UN: 3288
ADR: 6.1,II

**Specifications:**

Assay	>99%
Loss on drying	<0.5%
Calcium (Ca)	<0.005%
Chloride (Cl)	<0.005%
Sulfate (SO ₄)	<0.02%