

CL00.1151

ANALYT

Potassium hexacyanoferrate(II).3aq a.r.

Kaliumhexacyanoferrat(II).3aq z.A.
Potassio esacianoferrato(II).3aq p.a.
Kaliumhexacyanoferraat(II).3aq p.a.
Potassium hexacyanoferrate(II).3aq p.a.
Potassio hexacianoferrato(II).3aq p.a.

**For laboratory use, ISO
99+% K₄(Fe(CN)₆).3H₂O**

Mol weight: 422.39 g/mol
Density: 1,85 g/cm³
HS Nr: 28372000
CAS Nr: 14459-95-1
UN:
ADR:

Specifications:

Assay	>99%
Cadmium (Cd)	<0.0005%
Copper (Cu)	<0.002%
Sodium (Na)	<0.01%
Lead (Pb)	<0.002%
Carbonate (CO ₃)	<0.002%
Chloride (Cl)	<0.005%
Sulfate (SO ₄)	<0.005%
Non Solubles in Water	<0.005%