

Product Specifications

CL00.1166**Potassium nitrite v.p.**

Kaliumnitrit Rst.
Kaliumnitriet z.z.
Potassium nitrite t.p.
Potassio nitrato e.p.
Potasio nitrito t.p.

For laboratory use, FCC
97+%KNO₂

Mol weight: 85.11 g/mol
Density: 1,92 g/cm³
HS Nr: 28341000
CAS Nr: 7758-09-0
UN: 1488
ADR: 5.1,II

**Specifications:**

Assay	>97%
Arsenic (As)	<0.0001%
Iron (Fe)	<0.001%
Lead (Pb)	<0.0005%
Chloride (Cl)	<0.01%
Sulfate (SO ₄)	<0.02%