

Product Specifications

CL00.1435

ANALYT

Sodium molybdate a.r.

Natriummolybdaat.2aq p.a.
Sodio molibdato.2aq p.a.
Sodium molybdate.2aq p.a.
Natriummolybdat z.A.
Sodio molibdato.2aq p.a.

**For laboratory use, Ph. Eur., ACS
99.5+% Na₂MoO₄·2H₂O**

Mol weight: 241.95 g/mol
Density: 3.28 g/cm³
HS Nr: 28417000
CAS Nr: 10102-40-6
UN:
ADR:

Specifications:

Assay	>99.5%
Ammonium (NH ₄)	<0.001%
Iron (Fe)	<0.001%
Lead (Pb)	<0.001%
Chloride (Cl)	<0.001%
Phosphate (PO ₄)	<0.0005%
Sulfate (SO ₄)	<0.005%
pH	: 7.0 to 10.5 (5% solution @ 25°C)
Insoluble Matter (Non Solubles)	<0.005%
Heavy Metals as Lead (Pb)	<0.0005%