

Product Specifications

CL00.1457**Sodium dihydrogen phosphate.2aq v.p.**

Natriumdiwaterstofphosphat.2aq Rst.
Natriumdiwaterstoffosfaat.2aq z.z.
Sodium dihydrogène phosphate.2aq t.p.
Sodio diidrogeno fosfato.2aq e.p.
Sodio dihidrógeno fosfato.2aq t.p.

**For laboratory use, DAB, Ph. Eur., BP, USP, FCC
98-100.5% NaH₂PO₄·2H₂O**

Mol weight: 156.01 g/mol
Density: 1,92 g/cm³
HS Nr: 28352200
CAS Nr: 13472-35-0
UN:
ADR:

Specifications:

Assay	: 98-100.5%
Loss on drying	=21.5-24%
Arsenic (As)	<0.0002%
Calcium (Ca)	<0.002%
Copper (Cu)	<0.001%
Iron (Fe)	<0.0005%
Magnesium (Mg)	<0.002%
Lead (Pb)	<0.0005%
Zinc (Zn)	<0.001%
Chloride (Cl)	<0.005%
Fluoride (F)	<0.001%
Sulfate (SO ₄)	<0.02%