

Product Specifications

CL00.1464

di-Sodium hydrogen phosphate.2aq v.p.

di-Sodium hydrogène phosphate.2aq t.p.
di-Natriumwaterstoffosfaat.2aq z.z.
di-Sodio hidrógeno fosfato.2aq t.p.
di-Natriumwasserstoffphosphat.2aq Rst.
di-Sodio idrogeno fosfato.2aq e.p.

**For laboratory use, DAB, Ph. Eur., BP, Ph. Franç.
98.5-101% Na₂HPO₄·2H₂O**

Mol weight: 177.99 g/mol
Density: 2,10 g/cm³
HS Nr: 28352200
CAS Nr: 10028-24-7
UN:
ADR:

Specifications:

Assay	: 98.5-101%
Sodium dihydrogen phosphate	<1%
Loss on drying	=19.5-21%
Arsenic (As)	<0.0002%
Iron (Fe)	<0.001%
Heavy Metals as Lead (Pb)	<0.001%
Chloride (Cl)	<0.005%
Sulfate (SO ₄)	<0.01%