

Product Specifications

CL00.1470**Sodium nitrite v.p.**

Sodio nitrito e.p.
Sodio nitrito t.p.
Natriumnitriet z.z.
Natriumnitrit Rst.
Sodium nitrite t.p.

**For laboratory use, DAB8, B.P.C., USP, FCC
99-101% NaNO₂**

Mol weight: 69.00 g/mol
Density: 2,1 g/cm³
HS Nr: 28341000
CAS Nr: 7632-00-0
UN: 1500
ADR: 5.1 (6.1),III

**Specifications:**

Assay	: 99-101%
Loss on drying	<0.25%
Arsenic (As)	<0.0002%
Calcium (Ca)	<0.005%
Lead (Pb)	<0.0005%
Heavy Metals as Lead (Pb)	<0.001%
Chloride (Cl)	<0.005%
Sulfate (SO ₄)	<0.01%