

# Product Specifications

**CL00.1470**

## **Sodium nitrite v.p.**

Sodio nitrito e.p.  
Sodio nitrito t.p.  
Natriumnitriet z.z.  
Natriumnitrit Rst.  
Sodium nitrite t.p.

**For laboratory use, DAB8, B.P.C., USP, FCC  
99-101% NaNO<sub>2</sub>**

Mol weight: 69.00 g/mol  
Density: 2,1 g/cm<sup>3</sup>  
HS Nr: 28341000  
CAS Nr: 7632-00-0  
UN: 1500  
ADR: 5.1 (6.1),III



### **Specifications:**

Assay	: 99-101%
Loss on drying	<0.25%
Arsenic (As)	<0.0002%
Calcium (Ca)	<0.005%
Lead (Pb)	<0.0005%
Heavy Metals as Lead (Pb)	<0.001%
Chloride (Cl)	<0.005%
Sulfate (SO <sub>4</sub> )	<0.01%