

Product Specifications

CL00.1485

ANALYT

Sodium tungstate.2aq a.r.

Natriumwolframat.2aq z.A.
Sodium tungstate.2aq p.a.
Sodio tungstato.2aq p.a.
Sodio tungstato.2aq p.a.
Natriumwolframat.2aq p.a.

For laboratory use, ACS
99+% Na₂WO₄.2H₂O

Mol weight: 329.86 g/mol
Density: 3,23 g/cm³
HS Nr: 28418000
CAS Nr: 10213-10-2
UN:
ADR:

**Specifications:**

Assay	>99%
Free Alkali	<0.05%
Molybdenum (Mo)	<0.001%
Iron (Fe)	<0.001%
Heavy Metals as Lead (Pb)	<0.001%
Chloride (Cl)	<0.003%
Sulfate (SO ₄)	<0.01%
Nitrogen (N)	<0.001%