

Product Specifications

CL00.2020**Tetrachloroethylene v.p.**

Tetrachlorethylen Rst.
Tetrachloorethyleen z.z.
Tetracloroetileno t.p.
Tétrachloroéthylène t.p.
Tetracloroetileno e.p.

For laboratory use
99+% C2Cl4

Mol weight: 165.83 g/mol
Density: 1,623 g/ml
HS Nr: 29032300
CAS Nr: 127-18-4
UN: 1897
ADR: 6.1,III

**Specifications:**

Assay	>99%
Identity Infrared	: conforms
Acidity	< 0.0005meq/g
Non Volatiles	< 0.0005%
Water	< 0.01%