

Product Specifications

CL00.2601**Silver nitrate v.p.**

Silbernitrat Rst.
Argent nitrate t.p.
Zilvernitraat z.z.
Argento nitrate e.p.
Plata nitrate t.p.

**For laboratory use, DAB, Ph. Eur., BP, Ph. Franç., USP
99.8-100.5% AgNO₃**

Mol weight: 169.87 g/mol
Density: 4,35 g/cm³
HS Nr: 28432100
CAS Nr: 7761-88-8
UN: 1493
ADR: 5.1,II

**Specifications:**

Assay	: 99.8-100.5%
Copper (Cu)	<0.0005%
Iron (Fe)	<0.0005%
Lead (Pb)	<0.002%
Nitrite (NO ₂)	<0.1%
Sulfate (SO ₄)	<0.01%