

Product Specifications

CL00.2614

ANALYT

Silver nitrate a.r.

Argent nitrate p.a.
Zilvernitraat p.a.
Argento nitrate p.a.
Plata nitrate p.a.
Silberniträt z.A.
Prata nitrate p.a.

**For laboratory use, ACS, ISO, Ph. Eur.
99.99+% AgNO₃**

Mol weight: 169.87 g/mol
Density: 4,35 g/cm³
HS Nr: 28432100
CAS Nr: 7761-88-8
UN: 1493
ADR: 5.1,II



Specifications:

Assay	>99.99%
Cadmium (Cd)	<0.0001%
Copper (Cu)	<0.0002%
Iron (Fe)	<0.0002%
Manganese (Mn)	<0.0005%
Lead (Pb)	<0.001%
Zinc (Zn)	<0.0001%
Chloride (Cl)	<0.0001%
Nitrite (NO ₂)	<0.05%
Sulfate (SO ₄)	<0.002%
Clarity of the Solution	: Passes Test
Free Acid	: Passes Test
Substances not Precipitated by Hydrochloric Acid	<0.01%