

Product Specifications

CL00.2713

ANALYT

Aluminium(III) oxide, basic a.r.

Aluminium(III)oxide, basisch p.a.
Alluminio(III) ossido, di base p.a.
Aluminio(III) oxido, básico p.a.
Aluminium(III)oxid, Basich z.A.
Aluminium(III) oxide, basique p.a.

Mol weight: 101.96 g/mol
Density: 3,94 g/cm³
HS Nr: 28181011
CAS Nr: 1344-28-1
UN:
ADR:

**For laboratory use, activated for chromatography, Brockmann grade I
99+% Al₂O₃**

Specifications:

Assay	>99%
Particle size	: 40 µm <5%
pH	: 9.5 +/- 1.0
Silicate (SiO ₂)	<0.03%
Iron(III)oxide (Fe ₂ O ₃)	<0.03%
Sodium Oxide (Na ₂ O)	>0.4%
Chloride (Cl)	>0.002%