

Product Specifications

CL02.0221

CHROMALYT

Water, HPLC grade

Wasser, HPLC grade
Water, HPLC grade
Acqua, HPLC grade
Eau, HPLC grade
Agua, HPLC grade
Agua, HPLC grade

**For laboratory use, Reag. Ph. Eur., Reag. USP, HPLC
H2O - LF 1 μ S/cm - 0.4 μ m filtrated**

Mol weight: 18.016 g/mol
Density: 1,00 g/ml
HS Nr: 28539010
CAS Nr: 7732-18-5
UN:
ADR:

Specifications:

Non Volatiles	<0.0005%
Conductivity	<1 μ S/cm
T200nm >	97%
T210nm >	98%
T254nm >	99%
T300nm >	99%