

# Product Specifications

**CL02.0221****CHROMALYT****Water, HPLC grade**

Wasser, HPLC grade  
Water, HPLC grade  
Acqua, HPLC grade  
Eau, HPLC grade  
Agua, HPLC grade  
Agua, HPLC grade

**For laboratory use, Reag. Ph. Eur., Reag. USP, HPLC  
H2O - LF 1  $\mu$ S/cm - 0.4  $\mu$ m filtrated**

Mol weight: 18.016 g/mol  
Density: 1,00 g/ml  
HS Nr: 28539010  
CAS Nr: 7732-18-5  
UN:  
ADR:

**Specifications:**

Non Volatiles	<0.0005%
Conductivity	<1 $\mu$ S/cm
T200nm > 97%	
T210nm > 98%	
T254nm > 99%	
T300nm > 99%	