

# Product Specifications

**CL02.1126**

AUTOMATEX

## Potassium iodide 10% solution

Kaliumjodide 10% oplossing  
Potasio yoduro 10% solución  
Kaliumiodid 10% Lösung  
Potassio ioduro 10% soluzione  
Potassium iodure 10% solution

**TMAR - For laboratory use**  
**100 g KI / l H<sub>2</sub>O**

Mol weight: 166.01 g/mol  
Density: 1,07 g/ml  
HS Nr: 28276000  
CAS Nr: 7681-11-0  
UN:  
ADR:



### Specifications: