

Product Specifications

CL02.2714

Ammonium lactate 400 g/l solution

Ammonio lattato 400 g/l soluzione
Amonio lactato 400 g/l solución
Ammoniumlactat 400 g/l Lösung
Ammoniumlactaat 400 g/l oplossing
Ammonium lactate 400 g/l solution

For laboratory use
400 g C₃H₉NO₃ / l H₂O

Mol weight: 107.11 g/mol
Density: 1,103 g/ml
HS Nr: 29181100
CAS Nr: 515-98-0
UN:
ADR:

Specifications: