

Product Specifications

CL03.0204

HYDREX

Buffer pH 10.00

Tampão pH 10.00

Buffer pH 10.00

Tampon pH 10.00

Puffer pH 10.00

Tampón pH 10.00

Tampone pH 10.00

**For laboratory use, pH measurements, traceable to SRM from NIST
pH 10.00 ($\pm 0.02/20^\circ\text{C}$) $\text{H}_3\text{BO}_3\text{-KCl-NaOH}$**

**Traceability: Reference standards acc. NIST SRM 189c, 188, 185i,
186lg, 186llg, 187e**

Mol weight:

Density: 1,01 g/ml

HS Nr: 38229000

CAS Nr:

UN:

ADR:

Specifications:

pH	= 10.24 @ 0°C
pH	= 10.10 @ 10°C
pH	= 10.00 @ 20°C
pH	= 9.95 @ 25°C
pH	= 9.90 @ 30°C
pH	= 9.82 @ 40°C
pH	= 9.75 @ 50°C
pH	= 9.68 @ 60°C
pH	= 9.62 @ 70°C
pH	= 9.55 @ 80°C
pH	= 9.49 @ 90°C