

Product Specifications

CL03.0218

HYDREX

Buffer pH 9.00

Tampão pH 9.00
Buffer pH 9.00
Tampone pH 9.00
Tampón pH 9.00
Puffer pH 9.00
Tampon pH 9.00

**For laboratory use, pH measurements, traceable to SRM from NIST
pH 9.00 ($\pm 0.02/20^\circ\text{C}$) H₃BO₃-KCl-NaOH**

**Traceability: Reference standards acc. NIST SRM 189c, 188, 185i,
186lg, 186llg, 187e**

Mol weight:
Density: 1,01 g/ml
HS Nr: 38229000
CAS Nr:
UN:
ADR:

Specifications:

pH	= 9.24 @ 0°C
pH	= 9.11 @ 10°C
pH	= 9.00 @ 20°C
pH	= 8.97 @ 25°C
pH	= 8.93 @ 30°C
pH	= 8.86 @ 40°C
pH	= 8.80 @ 50°C
pH	= 8.75 @ 60°C
pH	= 8.71 @ 70°C
pH	= 8.67 @ 80°C
pH	= 8.64 @ 90°C