

Product Specifications

CL03.0706

HYDREX

Buffer pH 2.00 (colored)

Buffer pH 2.00 (gekleurd)
Tampão pH 2.00 (colorido)
Tampón pH 2.00 (de color)
Tampon pH 2.00 (coloré)
Tampone pH 2.00 (colorata)
Puffer pH 2.00 (farbcodiert)

**For laboratory use, pH measurements, traceable to SRM from NIST
pH 2.00 ($\pm 0.02/20^\circ\text{C}$) KCl-HCl (red)**

**Traceability: Reference standards acc. NIST SRM 189c, 188, 185i,
186lg, 186llg, 187e**

Mol weight:
Density: 1,00 g/ml
HS Nr: 38229000
CAS Nr:
UN:
ADR:

Specifications:

pH	= 1.99 @ 0°C
pH	= 1.99 @ 10°C
pH	= 2.00 @ 20°C
pH	= 2.00 @ 25°C
pH	= 2.00 @ 30°C
pH	= 2.00 @ 40°C
pH	= 2.00 @ 50°C
pH	= 2.00 @ 60°C
pH	= 2.00 @ 70°C
pH	= 2.00 @ 80°C
pH	= 2.00 @ 90°C