

**CL03.0710**

HYDREX

**Buffer pH 10.00 (colored)**

Tampão pH 10.00 (colori)  
Buffer pH 10.00 (gekleurd)  
Tampon pH 10.00 (coloré)  
Puffer pH 10.00 (farbcodiert)  
Tampón pH 10.00 (de color)  
Tampone pH 10.00 (colorata)

**For laboratory use, pH measurements, traceable to SRM from NIST  
pH 10.00 ( $\pm 0.02/20^{\circ}\text{C}$ ) H<sub>3</sub>BO<sub>3</sub>-KCl-NaOH (blue)**

**Traceability: Reference standards acc. NIST SRM 189c, 188, 185i,  
186lg, 186llg, 187e**

Mol weight:  
Density: 1,01 g/ml  
HS Nr: 38229000  
CAS Nr:  
UN:  
ADR:

**Specifications:**

pH	= 10.24 @ 0°C
pH	= 10.10 @ 10°C
pH	= 10.00 @ 20°C
pH	= 9.95 @ 25°C
pH	= 9.90 @ 30°C
pH	= 9.82 @ 40°C
pH	= 9.75 @ 50°C
pH	= 9.68 @ 60°C
pH	= 9.62 @ 70°C
pH	= 9.55 @ 80°C
pH	= 9.49 @ 90°C