

**CL03.0711**

HYDREX

**Buffer pH 9.00 (colored)**

Puffer pH 9.00 (farbcodiert)

Tampão pH 9.00 (colorido)

Buffer pH 9.00 (gekleurd)

Tampon pH 9.00 (coloré)

Tampón pH 9.00 (de color)

Tampone pH 9.00 (colorata)

**For laboratory use, pH measurements, traceable to SRM from NIST  
pH 9.00 ( $\pm 0.02/20^{\circ}\text{C}$ ) H<sub>3</sub>BO<sub>3</sub>-KCl-NaOH (green)****Traceability: Reference standards acc. NIST SRM 189c, 188, 185i,  
186lg, 186llg, 187e**

Mol weight:

Density: 1,00 g/ml

HS Nr: 38229000

CAS Nr:

UN:

ADR:

**Specifications:**

pH	= 9.24 @ 0°C
pH	= 9.11 @ 10°C
pH	= 9.00 @ 20°C
pH	= 8.97 @ 25°C
pH	= 8.93 @ 30°C
pH	= 8.86 @ 40°C
pH	= 8.80 @ 50°C
pH	= 8.75 @ 60°C
pH	= 8.71 @ 70°C
pH	= 8.67 @ 80°C
pH	= 8.64 @ 90°C