

Product Specifications

CL03.1112

HYDREX

Potassium nitrate 10% solution

Potasio nitrato 10% solución
Potassio nitrato 10% soluzione
Kaliumnitrat 10% Lösung
Kaliumnitraat 10% oplossing
Potassium nitrate 10% solution

Electrolite for pH electrodes - Filling solution
100 g KNO₃ / l H₂O

Mol weight: 101.11 g/mol
Density: 1,06 g/ml
HS Nr: 28342100
CAS Nr: 7757-79-1
UN:
ADR:

Specifications: