

Product Specifications

CL03.1124

HYDREX

Potassium chloride 0.2 mol/l (Conductivity standard)

Kaliumchloride 0.2 mol/l
Potassium chlorure 0.2 mol/l
Kaliumchlorid 0.2 mol/l
Potasio cloruro 0.2 mol/l
Potassio cloruro 0.2 mol/l

For laboratory use, LF measurements
LF Standard 22.4 mS/cm @ 20°C - 24.7 mS/cm @ 25°C (± 1%)
Traceability: Reference standards acc. NIST SRM 999b

Mol weight: 74.56 g/mol
Density: 1,01 g/ml
HS Nr: 28273985
CAS Nr: 7447-40-7
UN:
ADR:

Specifications: