

Product Specifications

CL05.0304

TITREX

Hydrochloric acid 0.01 mol/l

Acido cloridrico 0.01 mol/l
Acido clorhídrico 0.01 mol/l
Acide chlorhydrique 0.01 mol/l
Salzsäure 0.01 mol/l
Chloorwaterstofzuur 0.01 mol/l

For laboratory use, Titrations, ASTM D-4047
0.3646 g HCl / l H₂O = 0.01 N ($\pm 0.00002/20^{\circ}\text{C}$)
Traceability: Reference standards acc. NIST SRM 84I

Mol weight: 36.46 g/mol
Density: 1,00 g/ml
HS Nr: 28061000
CAS Nr: 7647-01-0
UN:
ADR:

Specifications: