

Product Specifications

CL05.0505

TITREX

EDTA diNa salt 0.05 mol/l

EDTA sal diNa 0.05 mol/l
EDTAdiNa zout 0.05 mol/l
EDTEdiNa salz 0.05 mol/l
EDTA sale diNa 0.05 mol/l
EDTA sel de diNa 0.05 mol/l

For laboratory use, Titrations

18.612 g EDTAdiNa.2H₂O / 1 H₂O = 0.05 M ($\pm 0.0002/20^{\circ}\text{C}$)

Traceability: Reference standards acc. NIST SRM 915b

Mol weight: 372.24 g/mol

Density: 1,01 g/ml

HS Nr: 29224985

CAS Nr: 6381-92-6

UN:

ADR:

Specifications: