

Product Specifications

CL05.1115

TITREX

Potassium thiocyanate 0.1 mol/l

Potassio tiocianato 0.1 mol/l
Kaliumthiocyanaat 0.1 mol/l
Potassium thiocyanate 0.1 mol/l
Kaliumthiocyanat 0.1 mol/l
Potasio tiocianato 0.1 mol/l

For laboratory use, Titrations
9.718 g KCNS / l H₂O = 0.1 N (±0.0002/20°C)
Traceability: Reference standards acc. NIST SRM 919b

Mol weight: 97.18 g/mol
Density: 1,01 g/ml
HS Nr: 28429080
CAS Nr: 333-20-0
UN:
ADR:



Specifications: