

# Product Specifications

**CL05.1115**

TITREX

## Potassium thiocyanate 0.1 mol/l

Potassio tiocianato 0.1 mol/l  
Kaliumthiocyanaat 0.1 mol/l  
Potassium thiocyanate 0.1 mol/l  
Kaliumthiocyanat 0.1 mol/l  
Potasio tiocianato 0.1 mol/l

**For laboratory use, Titrations**  
**9.718 g KCNS / l H<sub>2</sub>O = 0.1 N ( $\pm 0.0002/20^{\circ}\text{C}$ )**  
**Traceability: Reference standards acc. NIST SRM 919b**

Mol weight: 97.18 g/mol  
Density: 1,01 g/ml  
HS Nr: 28429080  
CAS Nr: 333-20-0  
UN:  
ADR:



### Specifications: