

Product Specifications

CL05.1116

TITREX

Copper(II) sulfate 0.01 mol/l

Cuivre(II) sulfate 0.01 mol/l
Rame(II) solfato 0.01 mol/l
Cobre(II) sulfato 0.01 mol/l
Kupfer(II)sulfat 0.01 mol/l
Koper(II)sulfaat 0.01 mol/l

For laboratory use, Titrations

2.497 g CuSO₄·5H₂O / l H₂O = 0.01 N (±0.00002/20°C)

Traceability: Reference standards acc. NIST SRM 3103a

Mol weight: 159.60 g/mol

Density: 1,01 g/ml

HS Nr: 28332500

CAS Nr: 7758-99-8

UN:

ADR:

Specifications: