

Product Specifications

CL05.1148

TITREX

Potassium hydroxide 1 mol/l in Ethanol

Kaliumhydroxide 1 mol/l in Ethanol
Potasio hidróxido 1 mol/l en Etanol
Potassium hydroxyde 1 mol/l dans l'Etanol
Kalilauge 1 mol/l in Ethanol
Potassio idrossido 1 mol/l in Alcole etilico

For laboratory use, Titrations (Solution may become brown on

standing !)

56.12 g KOH / l C₂H₅OH = 1 N (±0.005/20°C)

Traceability: Reference standards acc. NIST SRM 723e

Mol weight: -
Density: 0,86 g/ml
HS Nr: 38229000
CAS Nr:
UN: 2924
ADR: 3 (8),II

**Specifications:**