

Product Specifications

CL05.1406

TITREX

Sodium chloride 0.1 mol/l

Natriumchloride 0.1 mol/l
Sodio cloruro 0.1 mol/l
Sodium chlorure 0.1 mol/l
Sodio cloruro 0.1 mol/l
Natriumchlorid 0.1 mol/l

For laboratory use, Titrations
5.844 g NaCl / l H₂O = 0.1 N (±0.0002/20°C)
Traceability: Reference standards acc. NIST SRM 919b

Mol weight: 58.44 g/mol
Density: 1,00 g/ml
HS Nr: 25010099
CAS Nr: 7647-14-5
UN:
ADR:

Specifications: