

Product Specifications

CL05.1408

TITREX

Sodium hydroxide 0.01 mol/l

Natriumhydroxide 0.01 mol/l
Sodio hidróxido 0.01 mol/l
Natronlauge 0.01 mol/l
Sodio idrossido 0.01 mol/l
Sodium hydroxyde 0.01 mol/l

For laboratory use, TMAR - Titrations
0.4 g NaOH / l H₂O = 0.01 N ($\pm 0.00002/20^{\circ}\text{C}$)
Traceability: Reference standards acc. NIST SRM 723e

Mol weight: 40.00 g/mol
Density: 1,00 g/ml
HS Nr: 28151200
CAS Nr: 1310-73-2
UN: 1824
ADR: 8,III

Specifications: