

Product Specifications

CL05.1411

TITREX

Sodium hydroxide 0.1 mol/l

Sódio hidróxido 0.1 mol/l
Natriumhydroxide 0.1 mol/l
Sodium hydroxyde 0.1 mol/l
Sodio hidróxido 0.1 mol/l
Sodio idrossido 0.1 mol/l
Natronlaugé 0.1 mol/l

For laboratory use, Titrations, ASTM D-3712, D-4047

4 g NaOH / l H₂O = 0.1 N (±0.0002/20°C)

Traceability: Reference standards acc. NIST SRM 723e

Mol weight: 40.00 g/mol
Density: 1,01 g/ml
HS Nr: 28151200
CAS Nr: 1310-73-2
UN: 1824
ADR: 8,III



Specifications: