

Product Specifications

CL05.1417

TITREX

Sodium hydroxide 0.5 mol/l

Sodio idrossido 0.5 mol/l
Sodio hidróxido 0.5 mol/l
Sodium hydroxyde 0.5 mol/l
Natronlauge 0.5 mol/l
Natriumhydroxide 0.5 mol/l

For laboratory use, Titrations
20 g NaOH / l H₂O = 0.5 N (±0.001/20°C)
Traceability: Reference standards acc. NIST SRM 723e

Mol weight: 40.00 g/mol
Density: 1,02 g/ml
HS Nr: 28151200
CAS Nr: 1310-73-2
UN: 1824
ADR: 8,III



Specifications: