

Product Specifications

CL05.1424

TITREX

Sodium hydroxide 4 mol/l

Natriumhydroxide 4 mol/l
Sodio hidróxido 4 mol/l
Sodio idrossido 4 mol/l
Natronlauge 4 mol/l
Sodium hydroxyde 4 mol/l

For laboratory use, Titrations - TMAR Ammonium-Phenol Det. acc.

SERES

160 g NaOH / l H₂O = 4 N ($\pm 0.008/20^{\circ}\text{C}$)

Traceability: Reference standards acc. NIST SRM 723e

Mol weight: 40.00 g/mol

Density: 1,15 g/ml

HS Nr: 28151200

CAS Nr: 1310-73-2

UN: 1824

ADR: 8,II



Specifications: