

# Product Specifications

CL05.1912

TITREX

## Nitric acid 4 mol/l

Acido nítrico 4 mol/l  
Salpeterzuur 4 mol/l  
Salpetersäure 4 mol/l  
Acide nitrique 4 mol/l  
Acido nítrico 4 mol/l

**For laboratory use, Titrations**

**252.04 g HNO<sub>3</sub> / l H<sub>2</sub>O = 4 N (±0.008/20°C)**

**Traceability: Reference standards acc. NIST SRM 84I**

Mol weight: 63.01 g/mol

Density: 1,07 g/ml

HS Nr: 28080000

CAS Nr: 7697-37-2

UN: 2031

ADR: 8,II



### Specifications: