

# Product Specifications

**CL05.2603**

TITREX

## Silver nitrate 0.1 mol/l

Argent nitrate 0.1 mol/l  
Zilvernitraat 0.1 mol/l  
Silberniträt 0.1 mol/l  
Prata nitrato 0.1 mol/l  
Argento nitrato 0.1 mol/l  
Plata nitrato 0.1 mol/l

**For laboratory use, Titrations**

**16.987 g AgNO<sub>3</sub> / l H<sub>2</sub>O = 0.1 N (±0.0002/20°C)**

**Traceability: Reference standards acc. NIST SRM 3151**

Mol weight: 169.87 g/mol  
Density: 1,01 g/ml  
HS Nr: 28432100  
CAS Nr: 7761-88-8  
UN: 3082  
ADR: 9,III



### Specifications: