

Product Specifications

CL05.2610

TITREX

Sulfuric acid 0.01 mol/l

Schwefelsäure 0.01 mol/l
Acido solforico 0.01 mol/l
Acido sulfúrico 0.01 mol/l
Zwavelzuur 0.01 mol/l
Acide sulfurique 0.01 mol/l

For laboratory use, Titrations & TMAR TA-TAC Det.
0.981 g H₂SO₄ / l H₂O = 0.02 N (±0.00004/20°C)
Traceability: Reference standards acc. NIST SRM 84I

Mol weight: 98.08 g/mol
Density: 1,00 g/ml
HS Nr: 28070000
CAS Nr: 7664-93-9
UN:
ADR:

Specifications: