

Product Specifications

CL05.2612

TITREX

Sulfuric acid 0.05 mol/l

Acido solforico 0.05 mol/l
Acido sulfúrico 0.05 mol/l
Schwefelsäure 0.05 mol/l
Zwavelzuur 0.05 mol/l
Acide sulfurique 0.05 mol/l

For laboratory use, Titrations, ASTM D-3228
4.905 g H₂SO₄ / l H₂O = 0.1 N (±0.0002/20°C)
Traceability: Reference standards acc. NIST SRM 84I

Mol weight: 98.08 g/mol
Density: 1,01 g/ml
HS Nr: 28070000
CAS Nr: 7664-93-9
UN: 2796
ADR: 8,II



Specifications: