

# Product Specifications

CL06.0229

ANALYT

## Bipyridine-(2,2') a.r. (Dipyridyl)

Bipiridina-(2,2) p.a  
Bipiridina-(2,2) p.a  
Bipyridin-(2,2') z.A. (Dipyridyl)  
Bipyridine-(2,2') p.a. (Dipyridyl)  
Bipyridine-(2,2') p.a. (Dipyridyl)

For laboratory use, ACS, ISO, reagent for Fe and Mo Det.  
99.5+% C<sub>10</sub>H<sub>8</sub>N<sub>2</sub>

Mol weight: 158.18 g/mol  
Density: 1 g/cm<sup>3</sup>  
HS Nr: 29333999  
CAS Nr: 366-18-7  
UN: 2811  
ADR: 6.1,III



### Specifications:

Assay	>99.5%
Sulfated Ash	<0.2%
Non Solubles in Hydrochloric Acid	: passes test
Melting Point	> 70°C
Clarity of the Solution	: passes test