

Product Specifications

CL06.0608

INDEX

Phenolphthalein 5 g/l in ethanol/water (1/1) indicator solution

Phenolphthalein 5 g/l in Ethanol/Wasser (1/1) Indikator Lösung
Phénolphthaleïne 5 g/l dans l'éthanol/eau (1/1) indicateur solution
Fenolftaleína 5 g/l in etanolo/acqua (1/1) indicatore soluzione
Fenolftaleine 5 g/l in ethanol/water (1/1) indicator oplossing
Fenolftaleína 5 g/l en etanol/agua (1/1) indicador solución

TMAR TA-TAC Det. acc. AFNOR T.90.036
5 g C₂₀H₁₄O₄ / l C₂H₅OH/H₂O (1/1)

Mol weight: -
Density: 0,89 g/ml
HS Nr: 38229000
CAS Nr:
UN: 1993
ADR: 3,II



Specifications: