

Product Specifications

CL06.0609

INDEX

Phenolphthalein 1 g/l in ethanol indicator solution

Fenolftaleine 1 g/l in ethanol indicator oplossing
Phénolphtaléine 1 g/l dans l'éthanol indicateur solution
Fenolftaleina 1 g/l in etanolo indicatore soluzione
Fenolftaleína 1 g/l en etanol indicador solución
Phenolphthalein 1 g/l in Ethanol Indikator Lösung

AFNOR Det. in H₂O - T.I. (pH 8.2-9.8 - colourless to red-violet)
1 g C₂₀H₁₄O₄ / l C₂H₅OH

Mol weight: -
Density: 0,79 g/ml
HS Nr: 38229000
CAS Nr:
UN: 1993
ADR: 3,II



Specifications: